

Active walking and standing system

Thomy-Walker

Accessories list







Foot guides including fixation clip with rear twist lock, size 19 to 45



T-Handle for trunk stabilisation



Steering unit... with connecting parts for swivelling wheels without T-handle or front axle



Therapy table (transparent) without holder, table made of polycarbonate, with reinforced table top and aluminium rod



Lateral guide handlesCan be mounted on the hip arch, facing backward or forward



Mount for T-Handle or therapy table



Hand cuff

For positioning the hand on the T-Handle with Velcro fastener



Lumbar belt

for improved thoracic posture, neoprene, with four Velcro fasteners



Flat underarm support

Used for individuals with very weak upper body posture



Thorax- or pelvic pad

made of neoprene, to attach over the belt, available in different widths: 15, 20, 25, 30 cm



Underarm pads

lateral side support (1 pair is supplied with the basic model)



Buttock pad

for improved pressure distribution



Thigh fixation elastic (neoprene), for use in cases of internal thigh rotation



Push bar detachable with plug-in-fastener



Head or neck support complete with mounting parts



Gait assistance barWalking aid, complete with clamp, for step initiation



Knee extension lock to achieve better leg extension

Backstop suitable for standard wheels

